

ALLSHENG™ Mini centrifuge



The ALLSHENG™ Mini centrifuges are designed for quick spins of 1.5 ml/ 2.0 ml microcentrifuge tubes and PCR tubes at your lab workstation. The ALLSHENG™ Mini centrifuges are conveniently to operate. Simply close the lid to spin and open the lid to stop. Safety features help ensure the rotor will quickly decelerate to a stop when the lid is open.

Features:

- A. Super-Low Operating noise;**
- B. Small Size, Large Capacity;**
- C. Tool-Free Rotor Replacing;**
- D. Three Centrifugal Speed for Your Choice;**
- E. Decelerates Quickly and Safely**




Technical Data:

Model	Mini-6K	Mini-7K	Mini-10K
Power requirement	30W	30W	40W
Speed	6,000rpm/2000xg	7,000rpm/2300xg	10,000rpm/5000xg
Capacity	6x 1.5/2.0ml 16x0.2ml 0.2 ml 2x PCR Strips	6x 1.5/2.0ml 16x0.2ml 0.2 ml 2x PCR Strips	6x 1.5/2.0ml
Dimensions	175x148x118mm	175x148x118mm	175x148x118mm
Weight	0.7kg	0.7kg	0.7kg
Electrical	AC100-240V	AC100-240V	AC100-240V
Warranty	1Year		

Order information:

Code No.	Description
AS-08060-00-B/R/G	Mini-6K with BLUE/RED/GREEN Lid
AS-08070-00-B/R/G	Mini-7K with BLUE/RED/GREEN Lid
AS-08080-00-B/R/G	Mini-10K with BLUE/RED/GREEN Lid



TopLab Company (HK) Limited



(852) 9297-9506



info@hktoplalab.com



www.hktoplalab.com



Room 1608, F/16, GrandTech Centre, 8 On Ping Street, Shatin, NT., Hong Kong