



Mini-15K High Speed Mini Centrifuge



The Mini-15K High Speed Mini Centrifuge is your reliable laboratory assistant, which is compact, powerful and economical.

Humanized design, graceful shape and special processing of low noise design create a quiet, comfortable experimental environment for you.

Features:

- A. Extremely quiet operation**
- B. Easy-to-use display of time and speed**
- C. Autoclavable rotor (121°C, 20min)**
- D. Chemical resistant inner centrifuge lid**



Technical Data:

Model	Mini-15K
Speed	14,500 rpm (up to 14,000 xg)
Capacity	12 x 1.5 / 2.0ml micro-centrifuge tubes
Acceleration	15 seconds
Timer	20sec~99min
Dimensions (W x D x H)	270 x 230 x155 mm
Weight	5.0kg
Electrical	AC220V or AC120V, 50/60Hz
Warranty	1 Year

Order information:


Code No.	Description
AS-08110-00	Mini-15K, High speed Mini Centrifuge



TopLab Company (HK) Limited



(852) 9297-9506



info@hktoplab.com



www.hktoplab.com



Room 1608, F/16, GrandTech Centre, 8 On Ping Street, Shatin, NT., Hong Kong