





Before & After Centrifugation



During Centrifugation



MPC-P25 Mini-Plate Centrifuge is the first patented centrifuge for plates with an unique designed horizontal rotor in the world. All reactants are ensured to be in the bottom of the wells after centrifuging. It helps improve PCR yield or other tests.

Features:

- A. Unique design for the horizontal rotor centrifuge (Patented)**
- B. Two operations make the tests more convenient: Short-spin and timer setting**
- C. Quickly and smoothly braking from max. speed**
- D. Most types of PCR plates (skirted, semi, and non) and Elisa plates are acceptable**

Technical Data:

Model	MPC-P25
Speed	2,200 rpm/ 480xg
Capacity	2 x PCR Plates / Elisa Plates 8 x 0.2ml PCR tubes 8 x 0.2ml PCR strips
Dimensions (W x D x H)	290 x 360 x 140mm
Weight	3.6 kg
Electrical	AC220V or AC120V, 50/60Hz
Warranty	1 Year

Order information:

Code No.	Description
AS-08090-00	MPC-P25 Plate Centrifuge