

MTV-1 Mini Vortexer



With a 3mm orbit and fixed speed of 2,500rpm, the MTV-1 mini vortexer can vortex even up to 15 ml tubes. The MTV-1 has a small footprint of less than 74 x 91 mm, allowing it to fit on crowded benches.

MTV-1 can be powered by 4 pieces of 1.5V double A battery, avoiding your laboratories out of electricity. This is the most creative feature of this mini vortexer. Reliable and creative are the two best words to describe MTV-1 mini vortexer.



Features:

- A. Mini Size, 74 x 91 mm wide only**
- B. Powerful vortexing for up to 15 ml tubes**
- C. Powered by 1.5V batteries**
- D. Three different colors for choice**



Technical Data:

Speed Range	800 to 2,500 rpm
Orbit	3mm
Motor Activation	Rotary switch/ Touch pressure
Dimensions (WxDxH)	74x 91x 84 mm
Weight	0.35kg
Electrical	DC6V, (4pcs 1.5V double A battery)
Warranty	1 Year

Order information:

Code No.	Description
MTV-1-G	MTV-1 Mini Vortexer, Gold
MTV-1-B	MTV-1 Mini Vortexer, Blue
MTV-1-P	MTV-1 Mini Vortexer, Purple



TopLab Company (HK) Limited



(852) 9297-9506



info@hktoplab.com



www.hktoplab.com



Room 1608, F/16, GrandTech Centre, 8 On Ping Street, Shatin, NT., Hong Kong