

Power B1000 Power Supply



Specifications:

Output voltage	10 – 999 V in 1v steps
Output current	4 – 999 mA in 1 mA steps
Output accuracy	±1%
Output Power	200W max
Timing	1 – 999 min fully adjustable
Display	3-digital LED
Type of output	Constant voltage or Constant current
Pause/resume	Yes

Conditions:

Safety features	No load detection, sudden load change detection, overload / short circuit protection.
Dimension (LxWxH)	31x24x8 cm
Weight	1.5kg
Certificate	CE

Order information:

Code	Description
PP-1159	Power B1000 High Voltage Power Supply, 999V / 999mA / 200W



TopLab Company (HK) Limited



(852) 9297-9506



info@hktoplab.com



www.hktoplab.com



Room 1608, F/16, GrandTech Centre, 8 On Ping Street, Shatin, NT., Hong Kong