

# Power Bs Power Supply



## Specifications:

<b>Output voltage</b>	10 – 300 V in 1v steps
<b>Output current</b>	10 – 600 mA in 1 mA steps
<b>Output accuracy</b>	±1%
<b>Output Power</b>	100W max
<b>Timing</b>	1 – 999 min fully adjustable
<b>Display</b>	3-digital LED
<b>Type of output</b>	Constant voltage or Constant current
<b>Pause/resume</b>	Yes
<b>Warranty</b>	1 Year

## Conditions:

<b>Safety features</b>	No load detection, sudden load change detection, overload / short circuit protection.
<b>Dimension (L×W×H)</b>	25×19×8 cm
<b>Weight</b>	1.2kg
<b>Certificate</b>	CE

## Order information:

Code No.	Description
PP-1151	Power BS Power Supply, 300V / 600mA / 100W



TopLab Company (HK) Limited




(852) 9297-9506



info@hktoplab.com



www.hktoplab.com



Room 1608, F/16, GrandTech Centre, 8 On Ping Street, Shatin, NT., Hong Kong