

Homogenizer Bioprep-24

TOOLOO

Bioprep-24 is your reliable lab partner for sample preparation and cell lysis. Using this high-throughput homogenizer to extract DNA, RNA or proteins from animals, human and plant tissues, micro-organisms, spores, bones in less than 35 seconds. Whoever requires genomic DNA, RNA or recombinant protein as a starting material will benefit from this unit system.

Features:

- A. Efficient bead homogenization**
- B. 50 programmable memory settings and user optional condition settings**
- C. Disposable tubes ensure no cross-contamination**
- D. For many samples—from soft cells to bones**
- E. Large Capacity: 24*2ml tubes or 12*5ml tubes**



Order information:

Code No.	Description
AS-13010-01	Rotor for 24x2.0ml
AS-13010-02	Rotor for 6x5.0ml
AS-13010-03	2.0ml tubes
AS-13010-04	5.0ml tubes
AS-13010-05	Φ1.4mm Silica beads
AS-13010-06	Φ2.8mm Silica beads
AS-13010-07	Φ3mm Metal beads
AS-13010-08	Φ6mm Metal bead
AS-13010-09	Φ1mm glass beads
AS-13010-10	Φ3mm glass beads

Technical Data:

Model	Bioprep-24
Display	OLED
Speed Range (m/s)	4.00m/s~7.00m/s increment 0.05m/s
Capacity	6x 0.5ml/ 1.5ml/ 2ml or 2x 5ml tubes
Cycle per program	10
Noise Level	<70 db
Dimensions (W×D×H)	280×360×385 mm
Weight	25kg
Electrical	AC100~240V, 50/60Hz
Warranty	1 Year

TOOLOO

TopLab Company (HK) Limited



(852) 9297-9506



info@hktoplab.com



www.hktoplab.com



Room 1608, F/16, GrandTech Centre, 8 On Ping Street, Shatin, NT., Hong Kong