



TOOLOO

TRM-7S Digital Roller is designed for the gentle, thorough mixing of a variety of specimens, including blood samples, viscous substance and liquid-solid suspensions. Various containers, such as 15ml/ 50ml centrifuge tubes, test tubes, and roller bottles are acceptable.

Used in a variety of applications, prevention of blood coagulation and immune precipitation, etc. TRM-7S Digital Roller is an ideal choice for research, clinical laboratories.



Features:

- A. Completely Mixing: Rocking and Rolling Motions
- B. Gentle Rocking and Rolling
- C. Digital Speed and Time Setting
- D. LCD Display Clearly
- E. Small Footprint, but Large Capacity

Technical Data:

Model	TRM-7S
Speed	Adjustable, 10- 99 rpm (+/- 10%)
Tilt Angle	+/- 5°
Motion	Simultaneous rolling & rocking
Operating Temp. Range	+4 to +40° C

Dimensions (W x L x H)	270x 216 x 54 mm
Number of Rollers:	7
Weight (kg)	4.8
Electrical	80-230V,50-60Hz 50w
Warranty	1 Year

Order Information:

Code No.	Description
8011560	TRM-7S Digital Roller