



The Magic Programmable Mixer owns advanced programming function that far surpass most of traditional lab-use rotators.

The Magic Programmable Mixer can be selected to perform a wide variety of mixing actions: 1. 360° vertical rotating- for larger particles; 2. Swinging-for culture medium and larger particles medium; 3. Oscillating- for small particles.

In addition, the Magic Programmable Mixer includes an assortment of vertical tube holders. The Mini model is supplied with holders for 8 x 1.5/ 2.0ml tubes, 2 x 15ml and 2 x 50ml tubes; the Large model, 20 x 1.5/ 2.0ml tubes, 8 x 15ml tubes and 2 x 50ml tubes. For higher capacity applications, more holders can be purchased as alternatives. All kinds of tube holders can be assembled discretionarily in one rack.

Small footprint to fit on any crowded benches. The mixer can be easily moved to any place in your laboratory, cold rooms, incubators or biological hoods, etc. It is your reliable lab partner in mixing jobs.



Features:

- A. 10 different programs for rotation and mixing;
- B. 360° rotation and 0-360° discretionary mixing angle is available;
- C. Cold room and incubator compatible;
- D. Accept tubes from 1.5ml to 50ml;
- E. Tube holder assembled discretionarily;
- F. Speed and timer Adjustable



8011252-1



8011252-2



8011252-3



8011252-4



8011252-5

Technical Data:

Speed	Variable, 1 to 99 rpm
Capacity, Mini Model	8 x 1.5/ 2.0ml tubes, 2 x 15ml and 2 x 50ml tubes
Capacity, Large Model	20 x 1.5/ 2.0ml tubes, 8 x 15ml tubes and 2 x 50ml tubes
Operating Temp. Range	+4 to +65° C
Weight	1.6kg / 2.2kg
Adaptor	DC 12V, 1000 MA
Warranty	1 Year

Order Information:

Code No.	Model	Description
8011252	Mini	316mm length Magic Mixer with standard rack (8 x 1.5/ 2.0ml tubes, 2 x 15ml and 2 x 50ml tubes)
8011254	Large	420mm length Magic Mixer with standard rack (20 x 1.5/ 2.0ml tubes, 2 x 15ml and 2 x 50ml tubes)